

Advantages to Converting to a Bayonet Cuff



Old Luer Fitting



New Bayonet Fitting

CASMED standardized all of our Vital Signs Monitors to the bayonet fitting.

Here are some benefits for converting your facility to bayonet fittings:

- Bayonet fittings are more secure.
 - The bayonet fitting is a locking fitting. The male slip luer presents the possibility of coming undone during BP measurement.
- Bayonet fittings reduce the risk of cross-connection.
 - Male and female luer fittings can be accidentally connected to IV lines.
- Bayonet cuffs are always kept on the shelf at CASMED.
 - Results in shorter lead times.
- Bayonet cuffs are shipped as the standard for all monitor customers.
 - Reduces errors in shipping.
- Bayonet fittings are the most popular connections in NIBP on the market today. The monitor end is available in metal or plastic.